UNISOKU 2017

Newsletter

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COMPANY PROFILE

TRADE NAME	UNISOKU CO., LTD.
CEO	Shoji Suruga
FOUNDATION	NOV. 1974
LOCATION	2-4-3 Kasugano, Hirakata, Osaka, Japan
CAPITAL	50,000,000 JPY
BUSINESS	Manufacturing and sales of our own UHV LT SPM and kinetic spectroscopy system, research and development
BANKS	Bank of Kyoto, The Bank of Tokyo-Mitsubishi UFJ, The Shoko Chukin Bank
EMPLOYEES	44

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History

Nov. 1974	Establishment of Union Sokki Co., Ltd. Started the sales of Polarimeters and Light Scattering Spectrometers.
Oct. 1981	Changed the trade name to UNISOKU. Started the sales of Stopped Flow and Laser Photolysis Systems.
Oct. 1985	Increased capital to 5,000,000 JPY.
Apr. 1986	Completed and started the sales of Scanning Tunneling Microscopes for use in the atmosphere
Aug. 1989	Productized and started the sales of Ultra High Vacuum STM System.
Oct. 1990	Increased capital to 10,000,000 JPY.
Oct. 1991	Productized and started the sales of Ultra Low Temperature Ultra High Vacuum STM System.
Oct. 1997	Increased capital to 25,000,000 JPY.
Oct. 1998	Relocation of the head office due to business expansion.
Aug. 1999	Productized and started the sales of 2K High Magnetic Field Ultra Low Temperature Ultra High Vacuum STM System
Oct. 2002	Productized and started the sales of 400mK High Magnetic Field Ultra Low Temperature High Vacuum STM System.
Oct. 2002	Entered into Korean market
Oct. 2003	Entered into U. S. market
Aug. 2004	Entered into Chinese market
Feb. 2006	Received Nanotech 2006 "Evaluation and Measurement Award" for Scanning Probe Microscopes with UNISOKU's own Technologies.
Dec. 2006	Received the Excellence Award in "Manufacturing Business Award Ceremony" from Hirakata City.
Mar. 2007	Productized and started the sales of 4 Probe SPM System by JST Development of Systems and Technology for Advanced Measurement and Analysis
Apr. 2007	Awarded as "Person of merit for new technology development" from Osaka Prefecture for Ultra Low Temperature SPM.
Sep. 2010	UNISOKU becomes a group company of Tokyo Instruments Inc.
Sep. 2011	Establishment of UNISOKU-TII Instruments Co., Ltd. in Beijing, China, as an overseas affiliate.
Oct. 2011	Entered into global market (Europe, Asia, Australia, etc.) Increased capital to 50,000,000 JPY.
Jun. 2012	Completion of the first factory due to business expansion.
Dec. 2012	Productized and started the sales of Ultra Low Temperature High Magnetic Field STM System with Dilution Refrigerator.
Mar. 2014	Chosen as a "GLOBAL NICHE TOP 100" company from the Ministry of Economy, Trade and Industry of Japan
Jul. 2016	Completion of the second factory expanded 2.5 times in productivity, improving working environment

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Our duty is to bring about rich and convenient society through our customers' research results, and to contribute to development of scientific technologies

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UNISOKU and TII Group Are Growing

Globalization

UNISOKU STM (Scanning Probe Microscope) has more than 60% global share. We deliver our products throughout the world. We always have in mind improving UNISOKU's brand and popularity.

We aim to be a growing corporation, advancing in sales, technology, and capital tie-up, cutting edge with globalization and IoT.

Numbers of research and study reports using UNISOKU products are published on "Nature", "Science", etc.

Our original products are highly evaluated globally. UNISOKU was awarded as GNT-100 (Global Niche Top 100) from Ministry of Economy, Trade and Industry of Japan.

Our first factory was built in 2012, and second factory was completed in July 2016. Our productivity of probe microscopes expanded 2.5 times, Improvement of working environment.

Endeavors

Promoting "Mutual Satisfaction", we are building up satisfaction, mutual prosperity, win-win, business partnership, cooperation relationship with our customers, fellow traders and employees.

"Creation" valuing industry-academy cooperation. "Speed" accelerating practical application. Global partnership, cooperation.

By cooperation of Tokyo Instruments, advantage in optical measurements, and UNISOKU, surface analysis, we are challenging for world's first measuring technology combining probe microscopes and optics.

We grow up adapting to changes and quickly catching up with the times

Advancing into growing market, atomic size spectroscopic imaging

Basic approaches

Our No.1 priority is after sales support and product warranty

TII group put together the effort for our customers' problem solving

We improve quality with standardization and modulization. We also improve customer satisfaction by the reliability of products.

Being aware of IoT, we will boost in productivity, technologies, cost cuts, reduction of waste, to meet the deadlines definitely

UNISOKU's financial affairs are highly evaluated in "safety" "profitability" and "growth". We are ranked as good-standing company

Management Philosophy

- Contributing to development of Science
- Aiming No. 1 in nanotech measurement


Management Policy

Speed
Challenge
Positive Thinking

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
Prof. Q. K. Xue group (K. Chang et al.), Tsinghua Univ., discovered that atomic-thick SnTe thin film exhibits stable in-plane ferroelectricity by STM/STS measurements. The Curie temperature of this ferroelectricity is 270 K for 1-unit-cell (UC) thick SnTe, and 2- to 4- UC SnTe films are found to show robust polarization even at room temperature. STS measurements at the domain edge showed different electrical property dependent on the polarization direction, which may enable a wide range of application such as memory devices.

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
Recently, new 3D topological material, Weyl semimetal, is attracting much attentions. International research group lead by Prof. Z. Hasan, Princeton University (H. Zheng et al.), revealed topological surface states of Weyl semimetal NbP for the first time by using STM/STS measurements. They observed electronic standing waves of surface electrons in the real space (Fig. A) and quasi particle interference patterns in the reciprocal space (Fig. B). These results lead the unified understanding of Weyl semimetal and development of new nanodevices.

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
Prof. L. He group (L. J. Yin et al.) from Beijing Normal University succeeded in direct imaging of topological edge states* at the domain wall in graphene bilayer. From STS maps, these states were found to be quite robust even in the high magnetic fields (8T). Their result may raise hopes of graphene-based electronics with ultra-low dissipation.

*Topological edge states: 1D or 2D gapless states located at the edge or surface despite the bulk gap states

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
There has been no evidence of existence of Majorana Fermions* with which application to quantum computer is expected. The research group of Prof. S.C. Li, Prof. F.C. Zhang, Prof. J.F. Jia in Collaborative Innovation Center of Advanced Microstructures (H.H. Sun et al.) first observed strong spin dependence of tunneling effect (Spin Selective Andreev Reflection), which is predicted for Majorana Fermions, in the vortex of an artificial topological superconductor Bi₂Te₃ / NbSe₂ hetero-structure by using spin-polarized STM/STS.

*Majorana Fermions: a fermion whose antiparticle is itself which is theoretically predicted 80 years ago

Figures Courtesy of Dr. H.H. Sun and Prof. J.F. Jia


Thank you for using UNISOKU systems!

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International research group (Y.S. Fu et al.) led by Dr. T. Hanaguri, RIKEN, developed a new method to measure the magnetic moment of massless Dirac electrons at the surfaces of topological insulators by Scanning Tunneling Spectroscopy. Surprisingly, in two materials, Bi₂Se₃ and Sb₂Te₂Se, the magnetic moments of the Dirac electrons are very different both in amplitude and direction, whereas the velocities are the same. This result opens up a new way to control the spin of the topological surface state.

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
Dr. M. Allan, Leiden University, and co-researchers (I. Battisti, K.M. Bastiaans et al.), visualized the inhomogeneous melting process of frozen electrons in the iridates compound Sr_{2-x}LaxIrO₄ by using spectroscopic-imaging STM (Fig A shows data, Fig B shows a schematic). They discovered that the melting leads to a pseudo gap* and electronic order as in the cuprate high-temperature superconductors, establishing these phenomena as generic of lightly doped 2D Mott insulators.

*Pseudo-gap: unidentifiable gap structure which appears above superconducting transition temperature.

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
Prof. Q.K. Xue group in Tsinghua University, has been leading a research of high-temperature superconductivity (high-TC SC) over 65 K in a monolayer FeSe film on SrTiO₃. In a recent paper by H. Ding et al., they reported on the observation of high-TC SC and magnetic vortices in single-unit-cell FeSe films on anatase TiO₂(001) substrate by LT-STM, and clarified the essential roles of interface-induced electron-phonon coupling and charge transfer (CT) in boosting the high-TC SC. Moreover, they excluded the interfacial oxygen vacancies, and instead proposed the interfacial band alignment as primary source for CT between the substrate and FeSe films.


Figures Courtesy of Dr. C.L. Song and Prof. Q.K. Xue

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Cuprate superconductors show electronic bipartite – superconductivity and charge order. The group of Dr. T. Machida of Riken and Tokyo University obtained STS image, applying 11T magnetic field to one of cuprate superconductors, Bi₂Sr₂CaCu₂O₈+σ. As a result, they found that the charge order is enhanced by the introduction of the vortex which suppresses superconductivity at the high energy states above the pseudogap, meaning those two are in competitive relationship. This result will lead the breakthrough for the superconductivity mechanism investigation.

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
Electron in graphene has a novel property such as massless electrons and high carrier mobility, so the realization of superconductive state of graphene device has been expected. Prof. Takahashi from Tohoku University, Prof. Hasegawa, University of Tokyo and their collaborators (S. Ichinokura et al.), prepared Li- and Ca-intercalated bilayer graphene, C_6LiC_6 , C_6CaC_6 , and they confirmed that C_6LiC_6 doesn't show no superconductivity while the resistance of C_6CaC_6 starts to decrease at 4K and drops to zero at 2K by in situ 4-point-probe method combined with UHV-STM, meaning the achievement of the superconductive state.

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
Prof. W. F. Wang (X. Zhang et al.) from Peking University observed that simple V shape H3PH molecules were formed into motifs of Sierpinski Triangle (ST) or two dimensional crystals through hydrogen bond on Au(111) surface under low temperature (4.3K) and UHV (10^{-10} torr) condition. They found that the structure of the molecules such as ratio of ST was controllable by tuning the adsorbate density, chemistry of the surface and the intermolecular interactions. For example, by replacing Au base with Ag base which has stronger bonding with molecules, the ratio between 3-fold nodes and 4-fold nodes changed from 0% to 80%.

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
Silicene, 2D allotrope of silicon similar to graphene, is expected to be a silicon-based Dirac fermion material for the next generation nano-electronic device, but strong interaction with metallic substrate has been the key hurdle preventing its applications. Dr. Y. Du from University of Wollongong, Prof. J. Zhao from Dalian University of Technology and collaborators succeeded in preparing quasi-freestanding silicene in which top-layer is isolated from the substrate by intercalating the oxygen into bilayer silicene on Ag (111) substrate. The result of UHV-LT STM and Raman spectroscopy suggested that top-layer silicene exhibits 1x1 honeycomb structure which is different from silicene buffer layer.

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
Prof. Z.C. Dong, Prof. J.G. Hou and collaborators (Y. Zhang et al.) at University of Science and Technology of China measured the STM induced luminescence spectra from artificially constructed 2-to-4 phthalocyanine molecular oligomers, and succeeded in observing peak splitting caused by intermolecular dipole-dipole coupling. Furthermore, they performed STML mapping for each peak, and visualized the spatial distribution of the excitonic coupling for the phthalocyanine oligomers with sub-nanometer resolution in real space. This experimental approach can be a new tool to study intermolecular interaction at the level of individual molecules.

Figures Courtesy of Prof. Z.C. Dong, Prof. J.G. Hou


Thank you for using UNISOKU systems!

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Assoc. Prof. Otte group (R. Toskovic et al.), Delft University of Technology, constructed atomic Co chains of effective $S=1/2$ spins realizing the Heisenberg XXZ model in a transverse field by utilizing atom manipulation capability of STM. Inelastic tunneling spectroscopy (IETS) with varying magnetic field from 0 T to 9 T at 330 mK and Spin-Polarized STM for finely arranged 1~9 Co atoms reveals a series of ground state crossings which, in the thermodynamic limit, would constitute a quantum transition. 4-point-probe method combined with UHV-STM, meaning the achievement of the superconductive state.

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
Dr. H. Kim and Prof. Y. Hasegawa, Tokyo University, developed a method to measure the electrical conductance at different locations on individual atoms (PRL 114, 206801, 2015), and investigated the origin of a moiré pattern observed on a Pb overlayer. They measured the point contact conductance on sites equivalent in an atomic surface lattice but different in the moiré periodical structure, and found that the conductance does not depend on the contrast of the moiré pattern, and revealed that the moiré contrast observed by STM originates from the barrier height difference at small bias voltages.

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
Prof. J. Abe group (T. Yamaguchi et al.), Aoyama Gakuin University, succeeded in fast negative photochromism* which has been thought difficult to be realized. They synthesized binaphthyl-bridged phenoxyl-imidazolyl radical complexes, BN-PIC and BN-PIC2, which transform from the colored form to the colorless form by light irradiation. The mechanism was elucidated by using a nanosecond flash photolysis system and a cryostat.

*Photochromism: light induced reversible isomerization between two forms with different color.

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
Prof. S. Fukuzumi, Meijo University, and his collaborators (T. Tsudaka et al.) examined photodynamics of a series of 9-substituted 10-methylacridinium ions by three transient absorption techniques (pump-probe method, RIPT method* and nanosecond flash photolysis). They elucidated that intramolecular electron transfer (ET) from the donor moiety to the excited state of the acridinium ion moiety generated a long-lifetime triplet ET state, which has both high oxidizing and reducing activity (photosynthetic capability). 9-Methyl-10-methyl-acridinium ion exhibited the best performance in terms of the lifetime and quantum efficiency of the triplet ET state.

* RIPT method: Randomly - Interleaved - Pulse - Train method

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
2016 Physical Science award and this prize is called Chinese Nobel prize and prize is 1 million dollars.
 Formed single crystal thin film by MBE, discovered Quantum Anomalous Hall Effect (QAHE) and single layer ferrous (FeSe) high temperature superconductor phenomenon for the first time, using UNISOKU UHV STM system.

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
Apr.

駿河社長、中国（合肥、南京、上海、北京）を表敬訪問
 President Suruga made a courtesy visit on customers in China (Hofei, Nanjing, Shanghai, and Beijing)

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

What's New?


Low Temperature (LT)-UHV/Tip Enhanced Raman Spectroscopy USM1400 - LT TERS




Under LT UHV condition, sub-nanometer scale's Raman imaging with improved S/N ratio and sensitivity are realized.

- Seamless interface to TII's Nanofinder FLEX™
- Spatial resolution for Raman image : <10nm
- High enhancement factor by using Ag tip.
- Single molecular measurements are available


Picosecond Transient Absorption Spectrophotometer picoTAS

- Measures wide time range from 100 ps ~ s
- Completely Covers Gap Time Region of 1 ns ~ 20 ns
- Asynchronous Operation
- Broadband range from VIS ~ NIR
- Removes Fluorescence Signal
- Avoids Sample Damage
- Solution, Thin Film
- Compact Footprint

CROSS WORD

ACROSS →

- 2 One of units for pressure or (answer of 13)
- 4 TAS represents ___ absorption spectrometer
- 7 Femto-___-nano-micro-milli
- 9 The previous version of Windows 7
- 10 Raw seafood on rice
- 12 Binnig and Rohrer received ___ Prize for the design of STM
- 13 UHV represents ultra high ___

DOWN ↓

- 1 Unisoku is promoting ___ Satisfaction
- 2 Pen-___-apple-pen
- 3 Unisoku USM series has a capability of ___ polarized STM
- 5 The highest mountain in Japan is Mt. ___
- 6 Good - better - ___
- 8 Unisoku is in ___ prefecture, Japan
- 11 Carbon nanotube

for answers, take a look at www.unisoku.com

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Dassai is a Japanese sake getting popular recently.

It is the ultimate sake, "polishing" rice to the extreme.

"polishing" is a process of cleaning the rice after milling to take off proteins and make the taste clearer. The more the rice is polished, more cost and time are required. But this process is vital for making this significant sake, Dassai.

Dassai is awarded a gold prize by Monde Selection!

The news told that Prime Minister Abe gave Dassai to President Putin and President Obama as a present.

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UNISOKU Co., Ltd.